



msrB Antibody

Product Code	CSB-PA485132XA01E00
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	B7M1J2
Immunogen	Recombinant Escherichia coli O8 (strain IA11) msrB protein
Raised In	Rabbit
Species Reactivity	Escherichia coli O8 (strain IA11)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Escherichia coli O8 (strain IA11)
Lead Time	Made-to-order (14-16 weeks)
Target Names	msrB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.