



mscS Antibody

Product Code	CSB-PA843544XA01HUS
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	T0DVE4
Immunogen	Recombinant Helicobacter pylori (Campylobacter pylori) mscS protein
Raised In	Rabbit
Species Reactivity	Helicobacter pylori (Campylobacter pylori)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Helicobacter pylori (Campylobacter pylori)
Lead Time	Made-to-order (14-16 weeks)
Target Names	mscS
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.