



Metallothionein type 2b Antibody

Product Code	CSB-PA219915XA01FVW
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q19LA2
Immunogen	Recombinant Colocasia esculenta (Wild taro) (Arum esculentum) Metallothionein type 2b protein
Raised In	Rabbit
Species Reactivity	Colocasia esculenta (Wild taro) (Arum esculentum)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Colocasia esculenta (Wild taro) (Arum esculentum)
Lead Time	Made-to-order (14-16 weeks)
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.