



murl Antibody

Product Code	CSB-PA550421XA01PED
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q08783
Immunogen	Recombinant Pediococcus pentosaceus murl protein
Raised In	Rabbit
Species Reactivity	Pediococcus pentosaceus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Pediococcus pentosaceus
Lead Time	Made-to-order (14-16 weeks)
Target Names	murl
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.