



nanH Antibody

Product Code	CSB-PA341212XA01BDQ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P31206
Immunogen	Recombinant Bacteroides fragilis (strain YCH46) nanH protein
Raised In	Rabbit
Species Reactivity	Bacteroides fragilis (strain YCH46)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Bacteroides fragilis (strain YCH46)
Lead Time	Made-to-order (14-16 weeks)
Target Names	nanH
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.