



nhaP2 Antibody

Product Code	CSB-PA485903XA01EOR
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	B7LSJ3
Immunogen	Recombinant Escherichia fergusonii (strain ATCC 35469 / DSM 13698 / CDC 0568-73) nhaP2 protein
Raised In	Rabbit
Species Reactivity	Escherichia fergusonii (strain ATCC 35469 / DSM 13698 / CDC 0568-73)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Escherichia fergusonii (strain ATCC 35469 / DSM 13698 / CDC 0568-73)
Lead Time	Made-to-order (14-16 weeks)
Target Names	nhaP2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.