



nirS Antibody

Product Code	CSB-PA303275XA01PCA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P72181
Immunogen	Recombinant Paracoccus pantotrophus (Thiosphaera pantotropha) nirS protein
Raised In	Rabbit
Species Reactivity	Paracoccus pantotrophus (Thiosphaera pantotropha)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Paracoccus pantotrophus (Thiosphaera pantotropha)
Lead Time	Made-to-order (14-16 weeks)
Target Names	nirS
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.