



NEU2 Antibody

Product Code	CSB-PA734101XA01DXU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q64393
Immunogen	Recombinant Cricetulus griseus (Chinese hamster) (Cricetulus barabensis griseus) NEU2 protein
Raised In	Rabbit
Species Reactivity	Cricetulus griseus (Chinese hamster) (Cricetulus barabensis griseus)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Cricetulus griseus (Chinese hamster) (Cricetulus barabensis griseus)
Lead Time	Made-to-order (14-16 weeks)
Target Names	NEU2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.