



nhaS4 Antibody

Product Code	CSB-PA441228XA01SSQ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P72973
Immunogen	Recombinant Synechocystis sp. (strain PCC 6803 / Kazusa) nhaS4 protein
Raised In	Rabbit
Species Reactivity	Synechocystis sp. (strain PCC 6803 / Kazusa)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Synechocystis sp. (strain PCC 6803 / Kazusa)
Lead Time	Made-to-order (14-16 weeks)
Target Names	nhaS4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.