



nhhB Antibody

Product Code	CSB-PA325036XA01RDH
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P21220
Immunogen	Recombinant Rhodococcus rhodochrous nhhB protein
Raised In	Rabbit
Species Reactivity	Rhodococcus rhodochrous
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Rhodococcus rhodochrous
Lead Time	Made-to-order (14-16 weeks)
Target Names	nhhB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.