



ninG Antibody

Product Code	CSB-PA526164XA01EKZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O48427
Immunogen	Recombinant Enterobacteria phage H19B (Bacteriophage H19B) ninG protein
Raised In	Rabbit
Species Reactivity	Enterobacteria phage H19B (Bacteriophage H19B)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Enterobacteria phage H19B (Bacteriophage H19B)
Lead Time	Made-to-order (14-16 weeks)
Target Names	ninG
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.