



# EFER\_1155 Antibody

|                            |  |
|----------------------------|--|
| <b>Product Code</b>        | CSB-PA486731XA01EOR  |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                      |
| <b>Uniprot No.</b>         | B7LP27   |
| <b>Immunogen</b>           | Recombinant Escherichia fergusonii (strain ATCC 35469 / DSM 13698 / CDC 0568-73) EFER_1155 protein |
| <b>Raised In</b>           | Rabbit   |
| <b>Species Reactivity</b>  | Escherichia fergusonii (strain ATCC 35469 / DSM 13698 / CDC 0568-73)                               |
| <b>Tested Applications</b> | ELISA, WB (ensure identification of antigen)   |
| <b>Form</b>                | Liquid   |
| <b>Conjugate</b>           | Non-conjugated   |
| <b>Storage Buffer</b>      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, pH 7.4                   |
| <b>Purification Method</b> | Antigen Affinity Purified  |
| <b>Isotype</b>             | IgG  |
| <b>Clonality</b>           | Polyclonal   |
| <b>Product Type</b>        | Polyclonal Antibody  |
| <b>Immunogen Species</b>   | Escherichia fergusonii (strain ATCC 35469 / DSM 13698 / CDC 0568-73)                               |
| <b>Lead Time</b>           | Made-to-order (14-16 weeks)  |
| <b>Target Names</b>        | EFER_1155  |
| <b>Usage</b>               | For Research Use Only. Not for use in diagnostic or therapeutic procedures.                        |