



nhaP2 Antibody

Product Code	CSB-PA486170XA01EOP
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	B7MTW9
Immunogen	A synthesized peptide derived from Escherichia coli O81 (strain ED1a) nhaP2 protein
Raised In	Rabbit
Species Reactivity	Escherichia coli O81 (strain ED1a)
Tested Applications	ELISA, Dot blot (guaranteed specific to the immunizing peptide)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Escherichia coli O81 (strain ED1a)
Lead Time	Made-to-order (14-16 weeks)
Target Names	nhaP2
Product Modify	①. 100ug * KLH-conjugated peptide paired with 100ug * BSA-conjugated peptide (positive control); ②. 1ml*pre-immune serum (used as negative control) for free; ③. Rabbit polyclonal antibodies purified by Antigen Affinity;
Product Precision	1. CUSABIO guarantees the antibody purity > 95% confirmed by SDS-PAGE. 2. CUSABIO guarantees the antibody titer > 1: 64,000 confirmed by ELISA. 3. CUSABIO guarantees to offer Peptide-specific Dot blot. 4. CUSABIO guarantees 100% risk-free: No qualified antibody, no charge.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.