



CYP55A1 Antibody

Product Code	CSB-PA941412XA01FSH
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P23295
Immunogen	Recombinant Fusarium oxysporum (Fusarium vascular wilt) CYP55A1 protein
Raised In	Rabbit
Species Reactivity	Fusarium oxysporum (Fusarium vascular wilt)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Fusarium oxysporum (Fusarium vascular wilt)
Lead Time	Made-to-order (14-16 weeks)
Target Names	CYP55A1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.