



NTHL1 Antibody

Product Code	CSB-PA108707XA01CH
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A7M7B9
Immunogen	Recombinant Gallus gallus (Chicken) NTHL1 protein
Raised In	Rabbit
Species Reactivity	Gallus gallus (Chicken)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Gallus gallus (Chicken)
Lead Time	Made-to-order (14-16 weeks)
Gene Names	NTHL1