



ompC Antibody

Product Code	CSB-PA358176XA01SWW
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0A264
Immunogen	Recombinant Salmonella typhi ompC protein
Raised In	Rabbit
Species Reactivity	Salmonella typhi
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Salmonella typhi
Lead Time	Made-to-order (14-16 weeks)
Target Names	ompC
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.