



acyl Antibody

Product Code	CSB-PA899110XA01FWC
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q05053
Immunogen	Recombinant Pseudomonas sp. (strain V22) acyl protein
Raised In	Rabbit
Species Reactivity	Pseudomonas sp. (strain V22)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Pseudomonas sp. (strain V22)
Lead Time	Made-to-order (14-16 weeks)
Gene Names	acyl