



Ornithine aminotransferase Antibody

Product Code	CSB-PA016242XA01VCU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P31893
Immunogen	Recombinant Vigna aconitifolia (Moth bean) (Phaseolus aconitifolius) Ornithine aminotransferase protein
Raised In	Rabbit
Species Reactivity	Vigna aconitifolia (Moth bean) (Phaseolus aconitifolius)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Vigna aconitifolia (Moth bean) (Phaseolus aconitifolius)
Lead Time	Made-to-order (14-16 weeks)
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.