



phaB Antibody

Product Code	CSB-PA343466XA01PBZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P50204
Immunogen	Recombinant Paracoccus denitrificans phaB protein
Raised In	Rabbit
Species Reactivity	Paracoccus denitrificans
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Paracoccus denitrificans
Lead Time	Made-to-order (14-16 weeks)
Target Names	phaB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.