



phIB Antibody

Product Code	CSB-PA323066XA01STI
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P18954
Immunogen	Recombinant Serratia liquefaciens phIB protein
Raised In	Rabbit
Species Reactivity	Serratia liquefaciens
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Serratia liquefaciens
Lead Time	Made-to-order (14-16 weeks)
Target Names	phIB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.