



pnl Antibody

Product Code	CSB-PA335018XA01FWG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P24112
Immunogen	Recombinant Pectobacterium carotovorum (Erwinia carotovora) pnl protein
Raised In	Rabbit
Species Reactivity	Pectobacterium carotovorum (Erwinia carotovora)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Pectobacterium carotovorum (Erwinia carotovora)
Lead Time	Made-to-order (14-16 weeks)
Target Names	pnl
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.