



pheA Antibody

Product Code	CSB-PA333272XA01FGJ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P27603
Immunogen	Recombinant Pseudomonas stutzeri (Pseudomonas perfectomarina) pheA protein
Raised In	Rabbit
Species Reactivity	Pseudomonas stutzeri (Pseudomonas perfectomarina)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Pseudomonas stutzeri (Pseudomonas perfectomarina)
Lead Time	Made-to-order (14-16 weeks)
Target Names	pheA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.