



Pyridoxal phosphate homeostasis Antibody

Product Code	CSB-PA347292XA01VEU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P52055
Immunogen	Recombinant Vibrio alginolyticus Pyridoxal phosphate homeostasis protein
Raised In	Rabbit
Species Reactivity	Vibrio alginolyticus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Vibrio alginolyticus
Lead Time	Made-to-order (14-16 weeks)
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.