



PLSC Antibody

Product Code	CSB-PA675965XA01FWH
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q42870
Immunogen	Recombinant Limnanthes douglasii (Douglas' meadowfoam) PLSC protein
Raised In	Rabbit
Species Reactivity	Limnanthes douglasii (Douglas' meadowfoam)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Limnanthes douglasii (Douglas' meadowfoam)
Lead Time	Made-to-order (14-16 weeks)
Gene Names	PLSC