



pepX Antibody

Product Code	CSB-PA871745XA01LNG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9CE01
Immunogen	Recombinant Lactococcus lactis subsp. lactis (strain IL1403) (Streptococcus lactis) pepX protein
Raised In	Rabbit
Species Reactivity	Lactococcus lactis subsp. lactis (strain IL1403) (Streptococcus lactis)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Lactococcus lactis subsp. lactis (strain IL1403) (Streptococcus lactis)
Lead Time	Made-to-order (14-16 weeks)
Target Names	pepX
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.