



pfkA Antibody

Product Code	CSB-PA685401XA01FWI
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q59214
Immunogen	Recombinant Paenibacillus macquariensis (Bacillus macquariensis) pfkA protein
Raised In	Rabbit
Species Reactivity	Paenibacillus macquariensis (Bacillus macquariensis)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Paenibacillus macquariensis (Bacillus macquariensis)
Lead Time	Made-to-order (14-16 weeks)
Target Names	pfkA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.