



Pigment production hydroxylase Antibody

Product Code	CSB-PA333183XA01RDJ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P26698
Immunogen	Recombinant Rhodococcus sp. (strain ATCC 21145) Pigment production hydroxylase protein
Raised In	Rabbit
Species Reactivity	Rhodococcus sp. (strain ATCC 21145)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Rhodococcus sp. (strain ATCC 21145)
Lead Time	Made-to-order (14-16 weeks)
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.