



Genome polyprotein Antibody

Product Code	CSB-PA325529XA01PNN
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17767
Immunogen	Recombinant Plum pox potyvirus (strain Rankovic) (PPV) Genome polyprotein protein
Raised In	Rabbit
Species Reactivity	Plum pox potyvirus (strain Rankovic) (PPV)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Plum pox potyvirus (strain Rankovic) (PPV)
Lead Time	Made-to-order (14-16 weeks)
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.