



prl1 Antibody

Product Code	CSB-PA357641XA01OES
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09319
Immunogen	Recombinant Oreochromis mossambicus (Mozambique tilapia) (Tilapia mossambica) prl1 protein
Raised In	Rabbit
Species Reactivity	Oreochromis mossambicus (Mozambique tilapia) (Tilapia mossambica)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Oreochromis mossambicus (Mozambique tilapia) (Tilapia mossambica)
Lead Time	Made-to-order (14-16 weeks)
Target Names	prl1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.