



pri Antibody

Product Code	CSB-PA018724XA01DKC
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P48249
Immunogen	Recombinant Dicentrarchus labrax (European seabass) (Morone labrax) pri protein
Raised In	Rabbit
Species Reactivity	Dicentrarchus labrax (European seabass) (Morone labrax)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Dicentrarchus labrax (European seabass) (Morone labrax)
Lead Time	Made-to-order (14-16 weeks)
Gene Names	pri