



prl2 Antibody

Product Code	CSB-PA362735XA01FVK
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09584
Immunogen	Recombinant Oncorhynchus keta (Chum salmon) (Salmo keta) prl2 protein
Raised In	Rabbit
Species Reactivity	Oncorhynchus keta (Chum salmon) (Salmo keta)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Oncorhynchus keta (Chum salmon) (Salmo keta)
Lead Time	Made-to-order (14-16 weeks)
Target Names	prl2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.