



PSRP1 Antibody

Product Code	CSB-PA324869XA01FKI
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19954
Immunogen	Recombinant Spinacia oleracea (Spinach) PSRP1 protein
Raised In	Rabbit
Species Reactivity	Spinacia oleracea (Spinach)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Spinacia oleracea (Spinach)
Lead Time	Made-to-order (14-16 weeks)
Target Names	PSRP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.