



Recombinant Glucose-6-phosphate isomerase (pgi)

Product Code	CSB-MP009717KBE
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P77877
Product Type	Recombinant Protein
Immunogen Species	Klebsiella oxytoca
Purity	≥85% (SDS-PAGE)
Sequence	KHFSTTAPEK NLPVLLALIG IWYNNFFGAE TEAILPYDQY MHRFAAYFQQ GNMESNGKYV DRNGNAVYDQ TGPWIWGE PG TNGQHAFYQL IHQGTKMVPC DFIAPAITQN PLSDHHPKLL SNFFAQTEAL AFGKSREVVE QEYRDQ GKDP AALEHVVPFK VFEGNRP
Source	Mammalian cell
Target Names	pgi
Protein Names	Recommended name: Glucose-6-phosphate isomerase Short name= GPI EC= 5.3.1.9 Alternative name(s): Phosphoglucose isomerase Short name= PGI Phosphohexose isomerase Short name= PHI
Expression Region	1-167
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.