



ngl Antibody

Product Code	CSB-PA321901XA01EBK
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P21163
Immunogen	Recombinant Elizabethkingia miricola (Chryseobacterium miricola) ngl protein
Raised In	Rabbit
Species Reactivity	Elizabethkingia miricola (Chryseobacterium miricola)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Elizabethkingia miricola (Chryseobacterium miricola)
Lead Time	Made-to-order (14-16 weeks)
Target Names	ngl
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.