



ptsG Antibody

Product Code	CSB-PA681982XA01FLK
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q57071
Immunogen	Recombinant Staphylococcus carnosus (strain TM300) ptsG protein
Raised In	Rabbit
Species Reactivity	Staphylococcus carnosus (strain TM300)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Staphylococcus carnosus (strain TM300)
Lead Time	Made-to-order (14-16 weeks)
Target Names	ptsG
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.