



pstB1 Antibody

Product Code	CSB-PA341994XA01BRJ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P46342
Immunogen	Recombinant Bacillus subtilis (strain 168) pstB1 protein
Raised In	Rabbit
Species Reactivity	Bacillus subtilis (strain 168)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Bacillus subtilis (strain 168)
Lead Time	Made-to-order (14-16 weeks)
Target Names	pstB1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.