




putP Antibody

Product Code	CSB-PA703599XA01FKZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q53584
Immunogen	Synthesized peptide derived from internal of Staphylococcus aureus putP protein (1-15aa)
Raised In	Rabbit
Species Reactivity	Staphylococcus aureus
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Staphylococcus aureus
Target Names	putP

Image



Dot blot analysis of peptides corresponding to the immunogen was probed with putP Antibody at a dilution of 1:1000. Blocking/Dilution buffer: 4% nonfat dry milk/PBS. Secondary antibody: Goat polyclonal to rabbit IgG at 1/40000 dilution. Exposure time: 3min.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.