



pqqE Antibody

Product Code	CSB-PA516438XA01RJD
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O33506
Immunogen	Recombinant Rahnella aquatilis pqqE protein
Raised In	Rabbit
Species Reactivity	Rahnella aquatilis
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Rahnella aquatilis
Lead Time	Made-to-order (14-16 weeks)
Target Names	pqqE
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.