



prtD Antibody

Product Code	CSB-PA340556XA01FRO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P23596
Immunogen	Recombinant Dickeya chrysanthemi (Pectobacterium chrysanthemi) (Erwinia chrysanthemi) prtD protein
Raised In	Rabbit
Species Reactivity	Dickeya chrysanthemi (Pectobacterium chrysanthemi) (Erwinia chrysanthemi)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Dickeya chrysanthemi (Pectobacterium chrysanthemi) (Erwinia chrysanthemi)
Lead Time	Made-to-order (14-16 weeks)
Target Names	prtD
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.