



PROK Antibody

Product Code	CSB-PA361972XA01FWM
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P06873
Immunogen	Recombinant Parengyodontium album (Tritirachium album) PROK protein
Raised In	Rabbit
Species Reactivity	Parengyodontium album (Tritirachium album)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Parengyodontium album (Tritirachium album)
Lead Time	Made-to-order (14-16 weeks)
Target Names	PROK
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.