



pyrB Antibody

Product Code	CSB-PA309845XA01VCO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P96174
Immunogen	Recombinant Vibrio sp. (strain 2693) pyrB protein
Raised In	Rabbit
Species Reactivity	Vibrio sp. (strain 2693)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Vibrio sp. (strain 2693)
Lead Time	Made-to-order (14-16 weeks)
Target Names	pyrB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.