



URA5 Antibody

Product Code	CSB-PA325622XA01SIW
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P18904
Immunogen	Recombinant Sordaria macrospora URA5 protein
Raised In	Rabbit
Species Reactivity	Sordaria macrospora
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Sordaria macrospora
Lead Time	Made-to-order (14-16 weeks)
Target Names	URA5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.