



Putative peroxiredoxin in rubredoxin operon Antibody

Product Code	CSB-PA340522XA01CMA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P23161
Immunogen	Recombinant Clostridium pasteurianum Putative peroxiredoxin in rubredoxin operon protein
Raised In	Rabbit
Species Reactivity	Clostridium pasteurianum
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Clostridium pasteurianum
Lead Time	Made-to-order (14-16 weeks)
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.