



recA Antibody

Product Code	CSB-PA823790XA01BIP
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8VS69
Immunogen	Recombinant Borrelia hermsii recA protein
Raised In	Rabbit
Species Reactivity	Borrelia hermsii
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Borrelia hermsii
Lead Time	Made-to-order (14-16 weeks)
Target Names	recA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.