



5 Antibody

Product Code	CSB-PA890506XA01FRX
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9T1U3
Immunogen	Recombinant Acyrthosiphon pisum secondary endosymbiont phage 1 (Bacteriophage APSE-1) 5 protein
Raised In	Rabbit
Species Reactivity	Acyrthosiphon pisum secondary endosymbiont phage 1 (Bacteriophage APSE-1)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Acyrthosiphon pisum secondary endosymbiont phage 1 (Bacteriophage APSE-1)
Lead Time	Made-to-order (14-16 weeks)
Target Names	5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.