



RNA-directed RNA polymerase beta chain Antibody

Product Code	CSB-PA357727XA01FWN
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09674
Immunogen	Recombinant Escherichia phage BZ13 (Bacteriophage BZ13) RNA-directed RNA polymerase beta chain protein
Raised In	Rabbit
Species Reactivity	Escherichia phage BZ13 (Bacteriophage BZ13)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Escherichia phage BZ13 (Bacteriophage BZ13)
Lead Time	Made-to-order (14-16 weeks)
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.