



recX Antibody

Product Code	CSB-PA336363XA01FRQ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P37866
Immunogen	Recombinant Acidithiobacillus ferrooxidans (Thiobacillus ferrooxidans) recX protein
Raised In	Rabbit
Species Reactivity	Acidithiobacillus ferrooxidans (Thiobacillus ferrooxidans)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Acidithiobacillus ferrooxidans (Thiobacillus ferrooxidans)
Lead Time	Made-to-order (14-16 weeks)
Target Names	recX
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.