



CRO Antibody

Product Code	CSB-PA360939XA01FWQ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P03036
Immunogen	Recombinant Enterobacteria phage 434 (Bacteriophage 434) CRO protein
Raised In	Rabbit
Species Reactivity	Enterobacteria phage 434 (Bacteriophage 434)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Enterobacteria phage 434 (Bacteriophage 434)
Lead Time	Made-to-order (14-16 weeks)
Gene Names	CRO
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.