




rexB Antibody

Product Code	CSB-PA361145XA01FSF
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P03759
Immunogen	Recombinant Escherichia phage lambda (Bacteriophage lambda) rexB protein (27-46aa)
Raised In	Rabbit
Species Reactivity	Escherichia phage lambda (Bacteriophage lambda)
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Escherichia phage lambda (Bacteriophage lambda)
Target Names	rexB

Image



Western Blot

Positive WB detected in: recombinant protein

All lanes: rexB Antibody at 1:1000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 30 kDa

Observed band size: 30 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.