



SAP5 Antibody

Product Code	CSB-PA672570XA01SCG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q41389
Immunogen	Recombinant Saponaria officinalis (Common soapwort) (Lychnis saponaria) SAP5 protein
Raised In	Rabbit
Species Reactivity	Saponaria officinalis (Common soapwort) (Lychnis saponaria)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Saponaria officinalis (Common soapwort) (Lychnis saponaria)
Lead Time	Made-to-order (14-16 weeks)
Target Names	SAP5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.